



InBody520 Quick Installation Guide

Contents in the box

When opening the box, make sure all the following items are inside.

A. Product Units


- 1 InBody520 unit
- 2 Power Adapter (12V, 3.5A) (1EA)
- 3 Power Cord (AC 250V 10A, 70inch) (1EA)
- 4 User's manual CD
- 5 Quick installation guide
- 6 InBody520 Poster (1EA)(Optional)
- 7 Results Sheet Box (2EA, 500sheets / 1EA)(Optional)



3. Balance the InBody520

Check the horizontal level with the level indicator and the leveling screws.

The level indicator is located under the base frame cover and the leveling screws are located at the bottom of the base frame. If the equipment is not horizontally leveled, use the 5 leveling screws to level the equipment.



7. Self-calibration

- 1 When the power is turned on, there is a signal sound and as the logo comes up on the screen, self-calibration starts automatically.
- 2 Self-calibration will last about 2 minute.



- 3 When self-calibration finishes, the InBody520 becomes ready for the measurement.

When connecting peripherals(printers and other optional devices) to InBody520, turn on the power of peripherals before turning on the InBody520. When turning the power off, turn off the InBody520 before turning off the peripherals. This process will minimize the harm to the equipment caused by electric shock.

Do not touch the base frame when turning on the equipment and while it is warming up. During self-calibration, the loadcell inside the base frame is calibrated under the zero weight condition. Applying pressure or weight to the base frame during self-calibration will result in an inaccurate calibration.

Installation Instructions

A. Workplace Requirements


- (1) Location: Indoor only. Any outdoor area where the equipment is to be located should meet all the environmental requirements.
- (2) Operation environment: 10℃~40℃(50°F~104°F), 30%~80% RH, 500~1060hPa
- (3) Adapter : Power Input AC 100~240V, 50/60Hz, 1.2A
Power Output DC 12V, 3.5A

B. Attention for unpacking and assembling


Be sure to read below suggestions carefully before assembling.

- (1) Place the InBody520 on the horizontal ground.
- (2) Don't move the InBody520 by holding the hand electrode joints or LCD part.

Workplace environment




Right posture for moving




4. Put the results sheets in the printer

Please refer to the printer manual




The results sheets placement is determined by the printer brand. Please refer to your printer manual.

1. Unpack the box and remove the pads



2. Raise the upper part

- 1 Take the InBody520 out of the box and place it on the floor.
- 2 Draw the 'locking pin' located on the bottom right.
- 3 Raise the upper part.



5. Connect the printer and the adapter

- 1 Connect the USB cable into the USB HOST port. Either of two ports is available for use.


When using the printer, make sure to connect the InBody520 with the printer and turn on the printer before you turn on the InBody520. If you'd like to connect the printer when the InBody520 is turned on, you should turn off the InBody520, connect the printer, and turn it on again.

- 2 Connect the adapter cable into the power input port.

Use the adapter provided by Biospace to connect with the power port. When you connect the adapter cable, insert the adapter cable tightly into the power input port.

6. Turn on the power switch


- 3 The power switch is located on the back side of the equipment. Please turn on the equipment.




InBody520 Repacking

Be sure to turn off the power switch and unplug the power cable before repacking. Be careful to avoid severe physical shock, jarring or other damage while repacking, especially with regard to the hand module and foot electrodes.


- 1 Turn off the power.
- 2 Remove all cables connected to the InBody520.
- 3 Place the rectangular pad on the base frame.
- 4 Fold the upper part of the body, gently pulling the 'locking pin' located on the bottom right side.



- 5 Place the InBody520 into the box



- 6 Cover the top pad and seal the box with tape.




Error Messages

The InBody520 displays an error message to alert the operator of problems and to recommend corrective action.


A. Remove any objects on the footboard

When pressure or weight is applied to the base frame during warm-up, this message appears. Turn the InBody520 off and on again after removing the material on the base frame. Please do not apply pressure or weight to the base frame during the warm-up.




B. Incorrect Personal data

If inputted subject data are beyond acceptable ranges, this error message appears. When inputting subject data, do not exceed the acceptable ranges for Age and Height. Please refer to Chapter 2, Section 5: "Personal Profile" in the User's Manual.



C. Wipe your hands and feet with electrolyte tissue

If this message still displays on the LCD after checking the subject's posture during the second trial, this subject's soles and palms are too dry to be measured precisely. In this case, a subject should wipe his/her hands and feet with a wet tissue and try again.



InBody520 Setup Menu

When pressing the SETUP button, the following screen is displayed.

A. How to Modify Settings

- (1) Use the direction buttons (▲,▼) to select an item among Environment, Region, Printer, Result Sheet, Result Option, Others and Interface. Then, the selected item's subcategories are displayed on the screen.
- (2) Use the ENTER button or the direction button (▶) to move to the subcategories. Then, use the direction buttons (▲,▼) to move to the subcategory you want to modify and press the direction (▶).
- (3) Use the direction button (▲,▼) to modify the settings.
- (4) If you want to modify more items, press the EXIT/MODE button and use the direction button (▲,▼) to move and select the item you wish to modify. When you have no more to modify, press the EXIT/MODE button to move to another category: Environment, Region, Printer, Result Sheet, Result Option, Others and Interface. And then press the EXIT/MODE button.
- (5) When the screen asks whether or not you want to save the modified contents, press the ENTER or EXIT/MODE button to quit the Setup menu.

B. Setup

1) Environment: Shows the setup status of the InBody520.

- (1) Environment Summary
 - Language: Displays the language setup
 - Results Sheet: Displays the number of sheets that will be printed (0-2).
 - BMI Basis: Displays the selected standard rage for BMI.
(WHO: 18.5~25kg/m², Asian: 18.5~23kg/m², User: 18.5~27kg/m²)
 - Weight Measurement: MANUAL or AUTOMATIC will be displayed.

- (2) Memory Status
 - Total: Total capacity in use
 - Program: Capacity in use for program
 - Data: Capacity in use for result storage

- (3) Version Information: Displays the program version of InBody520.

2) Region

- (1) Set Data
 - Input the exact date you want to modify.
- (2) Set Time
 - Input the exact time you want to modify. It is 00(hour):00(minute):00(second) order.
- (3) Display Mode
 - Select a time display mode. (yy:year, mm:month, dd:day)
- (4) Unit
 - Select the unit for height and weight of personal information. (kg/cm, kg/in., lb./cm, lb./in.)
- (5) Language
 - Select the language you use. (English, Japanese, Korean, Spanish, Chinese, French, German, Italian, Russian)
- (6) Ethnic Background
 - Select the Ethnic Background of subject. (Asian, Caucasian, Hispanic, African, Others)


3) Printer

- (1) Printer: Select a printer you use.
- (2) Alignment
 - You can adjust the printing position. After adjusting, you can check whether the alignment has been set properly or not by printing a sample results sheet.
(Alignment range: X(left, right), Y(up, down) +50~-50)
- (3) Test Print
 - Check the printing position by printing out a sample results sheet. For this, move the selection box onto the Test Print icon using the direction buttons and press the ENTER button.

4) Result Sheet

- (1) Mode
 - Printed: This section allows you to print with Biospace's results sheet.
 - Built-in: This section allows you to print with plain white paper (A4 size).

- (2) Number of Result Sheet Printing
 - Select the number of sheets that will be printed (0-2).

 If it is set up as 'None', it will not print the results sheet.

5) Result Option

- (1) BMI Standard: Select the standard range for BMI.
 - 1 WHO: The standard range is 18.5~25 kg/m².
 - 2 Asian: The standard range is 18.5~23 kg/m².
 - 3 User: The standard range is 18.5~27kg/m².
- (2) Weight Control
 - 1 Enable: Print the result of Weight Control on the results sheet.
 - 2 Disable: No print the result of Weight Control.
- (3) Comprehensive Check
 - 1 Enable: Print the Comprehensive Check results (Nutritional Evaluation, Weight Management, Obesity Diagnosis, Body Balance, Body Strength) on the results sheet.
 - 2 Disable: No print the Comprehensive Check results.

6) Others

- (1) Measure Weight
 - 1 AUTO: When stepping on the base frame, weight is automatically measured and added to the personal information window.
 - 2 MANUAL: A subject directly enters his/her weight.
- (2) Adjust Weight

To perform this option use the up/down buttons (▲,▼) to adjust the weight. It is possible to calibrate the measured weight. If weight adjustment is needed due to clothes or accessories, select how much you want to default. (Calibration range: +5kg~-5kg, step: 0.1kg)

- (3) Adjust Volume
 - Adjust volume of the sound.
- (4) Sound Type
 - Select the sound type: Type1, Type2, Type3.
- (5) Initialize History
 - This is an option where you can erase all of the saved data from the InBody520. The history data on the InBody520 will be erased by pressing the ENTER button.
- (6) Gender Default
 - Select the gender automatically added to the personal information window. (Female, Male, Last Gender)

7) Interface

Setup the exterior device to connect to the InBody520 through a network.

- (1) Manual
 - Setting up the network is the same as setting up a network on the PC. (DNS, Netmask, Gateway, IP, Host IP)
- (2) Lookin' Body(PC, Data Management Program)
 - 1 Serial: Connecting the InBody520 to a PC with a Serial Cable to enable the use of Lookin'Body.
 - 2 USB: Connecting the InBody520 to a PC with a USB Cable to enable the use of Lookin'Body.
 - 3 Ethernet: Connecting the InBody520 to a PC with a Lan Cable to enable the use of Lookin'Body.
 - 4 Disable: Not using Lookin'Body.
- (3) Stadiometer
 - Select the application of the BSM330. (Enable/Disable)
- (4) Blood Pressure
 - Select the application of the TM2655, the blood pressure monitor. (Enable/ Disable)
- (5) DHCP (Dynamic Host Configuration Protocol)
 - DHCP offers IP address and dynamic Configuration of related data. When you use the network which uses DHCP server, get instructions from a system administrator.

C. Modification Example

- (1) In the case of measuring weight automatically.
 - Press the SETUP button.
 - Using the direction button (▲,▼), select the OTHERS menu.
 - Use the ENTER or direction button (▶) to move to MEASURE WEIGHT, and press the direction button (▶). Using the direction button (▲,▼), select AUTO.
 - Using the EXIT/MODE button, move to MEASURE WEIGHT and press the EXIT/MODE button twice to quit the SETUP menu and then press ENTER. Then the modified setup will be stored.
- (2) In the case of entering the weight directly.
 - Press the SETUP button.
 - Using the direction button (▲,▼), select the OTHERS menu.
 - Use the ENTER or direction button (▶) to move to MEASURE WEIGHT, and press the direction button (▶). Using the direction button (▲,▼), select MANUAL.
 - Using the EXIT/MODE button, move to MEASURE WEIGHT and press the EXIT/MODE button twice to quit the SETUP menu and then press ENTER. Then the modified setup will be stored.

InBody520 Quick Setup



Users can easily modify frequently changed functions in the initial screen.

A. How to modify using Quick Setup

- (1) Press the EXIT/MODE button on the initial screen.
- (2) Choose an item you want to modify.
 - Use the right and left direction button (◀▶) to move between items.
 - Use the up and down direction button (▲▼) to modify setup values.
- (3) Modified setup values are automatically stored. After finishing modification, press EXIT/MODE or the ENTER button to exit the Quick Setup.



B. Quick Setup

- (1) Weight Adjust
 - It is possible to calibrate the measured weight. If weight adjustment is needed due to clothes or accessories, select how much you want to deduct.
(Calibration range : +5kg ~ -5kg, step : 0.1kg)
- (2) Unit
 - Select the unit for height and weight of personal information. (kg/cm, kg/in., lb./cm, lb./in.)
- (3) Result Sheet Mode
 - Used to select the type of result sheet.
 Printed : Select when using the printed result sheet provided by Biospace.
 Built-in : Select when using plain white paper (A4 size).
- (4) Printing number of result sheets
 - Decide on the number of result sheets that will be printed.
- (5) Volume
 - Used to control the sound volume.
- (6) Date / Time
 - Used to set up the current time.


InBody520 Database





InBody520 has a data saving function in itself. When you input your I.D. the measured data is saved automatically. You can search, delete, back up and restore the data at the DATABASE screen and copy them using USB storage. The copied data can be checked at the excel file in computer.


A. DB MENU (Backup, Restore, Copy)

InBody520 provides whole DB data Backup, Restore and Copy functions.

Press the DATABASE button and the DATABASE screen appears.



- (1) Backup
 - This is for the backup of whole data saved in DB. Connect the USB to InBody520 and press the number '1' button on the keypad then it starts to backup.
 This Backup function is only for restoration using USB storage when needed. The backup data can not be identified in your computer. If you want to check the measurement data in your computer, use 'Copy' function.
- (2) Restore
 - This is for the restoration of the back up data to InBody520 again. Connect the USB storage that has backup data to InBody520 and press the number '2' button on the keypad then it starts the restoration.
 Note that stored data in the InBody520 is automatically deleted when restoring new data from USB.
 There is no limit to individual data storage. It is possible to store approximately 10,000 results in total.
- (3) Copy
 - This is for the copy of whole data in DB. You can check the copied data in Lookin'Body. If Lookin'Body is not installed, you can check the data in Excel file.
 Contact Biospace or an authorized distributor for information about USB storage devices compatible with InBody520

 How to check the copied data in computer

1) Check in excel file

Connect the USB that has backup data to the computer. You can find INBODY folder is created. In 'InBody' folder, Lookin'Body folder and a file which is '.csv' (db_backup.csv) are found. Open the db_backup.csv file and you can check the data in Excel file.


2) Check in Lookin'Body

In case of using Lookin'Body, use 'USB_Convert' (Start > Program > Lookin'Body3.0 > USB_Convert). Open USB_Convert and click 'Move to Lookin'Body' button then the data in USB storage moves to Lookin'Body automatically. After completion of transmission, you can search and check the data on Lookin'Body screen.



B. Data search

Press the DATABASE button and the DATABASE screen appears. Basically, the latest measured data for a month is shown.



There are two ways of data searching: ID search and date search

- (1) I.D. Search
 - 1 12+ENTER: to find all I.D.'s with 12 in the first 2digit.
 - 2 D+ENTER: to find all I.D.'s with D in the first digit.
 - 3 ENTER: to list all stored data.
- (2) Date Search
 - 1 Input the first day of the period you want to search (From).
 - 2 Input the last day of the period you want to search (To).
 - 3 Press ENTER button.
If you want to search the special date, input the same date in the both of From and To column.

C. Searched data management (Reprint, Copy, Delete)


You can reprint, delete and copy the results searched by I.D. or date. This is the sample screen of the searched data.



- (1) Reprint
 - Locate the selection box on the data you want to print using the direction button (▲,▼) and press number '1' button on the keypad to print the data.
- (2) Copy
 - Locate the selection box on the data you want to copy and press the ENTER button. If you want to copy more than 1, move the selection box using the direction button (▲,▼) and press ENTER. Select all data you want to copy and press the number '2' button, then the data will be copied into a USB storage device.



- (3) Delete
 - Locate the selection box on the data you want to delete and press the ENTER button. If you want to delete more than 1, move the selection box using the direction button (▲,▼) and press ENTER. Select all data you want to delete and press the number '3' button,


 Once deleted, data cannot be restored.

InBody520 Maintenance

- (1) After usage wipe electrodes with a wet tissue (also known as wet naps, surface wipes, anti-bacterial tissues, etc.). Wet tissues can be purchased at any local stores.
- (2) Handle the arms with care, Never apply excessive stress near the hand electrode joint. The damage caused by misuse may affect the function of the internal cable and electric board.
- (3) Do not place anything on the base frame nor apply any pressure onto it when the InBody520 is not in use.
- (4) When the InBody520 is not in use for a period longer than a day, unplug the adapter.
- (5) Do not move or relocate the InBody520 while the power is on.
- (6) Do not drop food or drinks on the equipment. They may affect the electrical parts in the equipment or cause damage.
- (7) Once a week, wipe the exterior sides of InBody520 with a dry towel. In particular, clean the LCD monitor, softly with care not to scratch the surface.
- (8) Follow local governing ordinances and recycling plans regarding disposal or recycling of device components.

Limited Warranties

- (1) Biospace guarantees that the product has been approved with the qualified test procedure under the severe condition.
- (2) The one year factory warranty begins on the date of purchase.
- (3) During the one - year warranty period, Biospace remedies any original defect in material or workmanship.
- (4) The following defects or malfunctions will not be covered under the one year warranty :
 - Any defect caused by user's fault.
 - Any defect or damage caused by not following the instructions described in the user's manual.
 - Any defect or damage caused by natural disasters (storm, flood, earthquake, etc.).
 - Any defect or damage caused by disassembly of InBody520 or by modifying internal parts or program by unauthorized person.
- (5) An extended warranty may be purchased by executing an extended warranty contract after the initial warranty period.
- (6) No return is allowed after a product is opened. The proof of purchase must be accompanied before requesting service.
- (7) Please contact an authorized service for any service calls.

 Contact BIOSPACE for details regarding the limited warranties.

Biospace

©1996-2007 Biospace Co., Ltd. All rights reserved.

No part of this Quick installation guide shall be reproduced, stored in any retrieval system, or transmitted by any means (electronic, mechanical, photocopying, recording, or otherwise) without written permission from Biospace. No patent liability is assumed with respect to the use of the information contained herein. This installation guide may contain misprints, which can be modified without prior notice to the readers. Biospace is not liable for any damage caused by a failure to meet the requirements in this installation guide.

Biospace Co. Ltd.

Tel : 82-2-501-3939 Fax : 82-2-501-3978

Homepage : <http://www.biospace.co.kr>

E-mail : biospace@biospace.co.kr

Biospace, the Biospace logo, InBody and Lookin' Body are registered trademarks of Biospace Co., Ltd. The names of the company and products in this installation guide other than those of Biospace Co., Ltd. are the trademarks of the company. Stating the products of other company is strictly for the purpose of providing information, not to guarantee or recommend the products. Biospace is not responsible for the performance or the use of these products.

Biospace reserves the right to modify the dimensions or exterior of InBody520 to improve the quality of the products, without consent of the customer.

InBody520